

VereFoodborne[™] Detection Kit Instructions for Use

For Research Use Only



VFBN-RA50



50

Store at -25°C to -15°C (frozen components) Store at 15°C to 25°C (ambient components)

83 Science Park Drive #04-02 Singapore (118258)

November 2018

IFU-FBN-RA01-1009



Contents

Purpos	е	2
Summa	ary and Explanation	2
Princip	e of the Procedure	3
Kit Con	tent	4
Storage	e Condition	4
Materia	Is Required but Not Provided	6
Warnin	gs and Precautions	7
Quality	Control	7
Workflo)W	8
Specim	en Collection, Handling and Storage	9
Brotoc	Preparation	9
Imports	ant notes before starting	10 10
1	Dilute sl. PCR Control	
2.	Prepare PCR Reaction Mix	
3.	Load PCR Reaction Mix into Chip	12
4.	Seal Chip for PCR	13
5.	Run Chip (PCR)	15
6.	Prepare Hybridization Mix	19
7.	Load Hybridization Mix into Chip	19
8.	Seal Chip for Hybridization	21
9.	Run Chip (Hybridization)	23
10.	Wash Chip	23
11.	Detection	25
Assay	Control	
Interpre	etation of Results	30
Trouble	eshooting Guide	31
Limitati	ons of the Test	
Dispos	al	
Contoo	cal Assistance	
Undors	l tanding the Symbols	ວວ ຈະ
Produc	t Use Limitations Warranty Disclaimer	
Notice	to Purchaser	38



Purpose

VereFoodborne[™] Detection Kit is a multiplex PCR/microarray-based research use only test. This nucleic acid-based test is intended for qualitative detection and differentiation of major food pathogens after nucleic acid extraction from a single source sample. One of the following organisms can be identified:

Vibrio cholera Vibrio parahaemolyticus Staphylococcus aureus Listeria spp. Bacillus spp. Clostridium perfringens Campylobacter jejuni/ Campylobacter lari/Campylobacter coli Shiga toxin-producing Escherichia Coli (STEC) Shigella other spp. Salmonella spp. Cronobacter sakazakii Shiga toxin genes (stx1A, stx2A)

This test is for use in conjunction with the VerePLEX™ Biosystem.

Summary and Explanation

Foodborne diseases are responsible for a wide range of illnesses and complications ranging from diarrheal diseases to various forms of cancer. Serious outbreaks of foodborne diseases have been documented in both developed and developing countries in the past decade, and the rate of illnesses are increasing significantly owing to the key factors of increased consumption of minimally processed food, the globalization of the food supply and the mass production and distribution of ready-to-eat food. In the United States alone, the Centers for Disease Control and Prevention (CDC) estimates that there are approximately 48 million cases of foodborne diseases, resulting in 128,000 hospitalizations and 3,000 deaths each year¹.

VereFoodborne[™] Detection Kit uses the lab-on-Chip (LOC) platform for simultaneous and qualitative detection and identification of major food pathogens based on selected target genes.

¹ Centers for Disease Control and Prevention (CDC). Estimates of Foodborne Illness in the United States. http://www.cdc.gov/foodborneburden/



Principle of the Procedure

VereFoodborne[™] Detection Kit includes a single-use disposable VereChip[™] as well as reagents and consumables necessary for nucleic acid amplification and DNA hybridization. PCR primers and probes are designed to target the different genes of target organisms as well as internal controls. This test is for use in conjunction with the VerePLEX[™] Biosystem.

The VerePLEX[™] Biosystem consists of the following components:

- Temperature Control System (TCS), which is a system that thermally drives the VereChip[™]. It consists of five Temperature Control Modules (TCMs) that allow five independent temperature programs to be run independently
- Optical Reader (OR), which detects and analyzes the microarray fluorescence in the range between 670 and 730 nm
- Touch Monitor (TOM), which connects to the TCS for display and input via touch
- Barcode Reader
- VerePLEX[™] Biosystem Software

The system requires a VereChip[™], on which a miniaturized reactor for PCR amplification, and a DNA microarray are integrated.

The DNA microarray consists of microscopic spots of DNA oligonucleotides called probes that are fixed on a silicon substrate. These probes are short sections of target genes which hybridizes to a complementary strand of amplified product. The amplified product is labeled with a fluorophore dye and probe-product hybridization is detected as a fluorescent signals and captured by CCD camera in the OR.

The probes for targets are spotted in duplicates in a 6 rows x 21 columns layout. If target pathogen is in the sample and captured, respective probes will light up in a particular pattern on the microarray. For a sample to be positive for a particular gene or pathogen, certain criteria must be satisfied. The criteria are written in "*Diagnostic Rule files*" provided in the BioApplication. The VerePLEXTM Biosystem Software relies on this set of rules for pattern interpretation.



Kit Content

Kit Catalog no.	VFBN-RA50	
Tests		50
Ambient Components (15°C to 25°C)	Part Number	Quantity
Box 1 of 2		
VereFoodborne™ Chips	FBN-30	2 boxes
PCR Clamps	CPR-01	1 pack
IN Clamps	CIN-02	2 packs
HYB Clamps	CHY-01	1 pack
Wash Buffer Concentrate	VCP-WBC03/ VCP-WBC03	1 bottle
Frozen Components (-25°C to -15°C)	Part Number	Quantity
Box 2 of 2		
VereFoodborne™ Primer Mix A	FBN-PMA04	2 tubes
VereFoodborne™ Primer Mix B	FBN-PMB03	2 tubes
VCP Hyb Buffer	VCP-15B01	2 tubes
Hyb Probe Concentrate	VCP-HPC01	2 tubes
sL PCR Control	VCP-SPC01	2 tubes



Storage Condition

- Store all Frozen Components at -25°C to -15°C upon receipt.
- VereFoodborne[™] Primer Mix A, VereFoodborne[™] Primer Mix B and sL PCR Control are stable for up to 2 months after opening. Keep the number of freezethaw cycles to <6 cycles.
- Store diluted sL PCR Control (5x10² copies/µL) at -25°C to -15°C for up to 1 month. Keep the number of freeze thaw cycles to < 3 cycles.
- Keep the Primer Mix and Hyb Probe Concentrate away from light until ready to use.

NOTE: Repeated thawing and freezing may reduce the sensitivity of the assay.

- Store reconstituted Microarray Hybridization Mix at 2°C to 8°C for up to 2 months.
- Store all Ambient Components at room temperature (15°C to 25°C).
- After opening the aluminum packaging, store the Chips in its original box at room temperature for up to 2 months.

NOTE: Keep Chips away from light.

- Precipitation or crystallization may occur in Wash Buffer Concentrate if the storage temperature is low. Should this occur, please refer to important note before Chip washing step (page 23).
- Once opened, Wash Buffer Concentrate can be kept at room temperature for up to 1 year or until expiration date, whichever comes first.
- Store diluted Wash Buffer at room temperature for up to 2 months. Precipitation or crystallization may occur in diluted Wash Buffer if the storage temperature is low. Should this occur, please refer to important note before Chip washing step (page 24).
- If left unopened, all reagents are stable until the expiration date indicated on the respective labels.



Materials Required but Not Provided

- Reagents a)
 - DNeasy *mericon* Food Kit² or equivalent •
 - QuantiTect Multiplex PCR NoRox Kit³ or • equivalent
 - PCR Grade Water
 - Distilled/Reverse Osmosis (RO)/Ultrapure Water
- Consumables b)
 - Personal protective equipment •
 - Sterile beveled pipette tips⁴ for Chip loading •
 - Sterile filter pipette tips •
 - 1.5 mL microcentrifuge tubes
 - 0.2 mL PCR tubes
 - Sterile 50 mL centrifuge tube (non-skirted)
 - 1 L Bottle
 - **Decontamination products** •



- Centrifuge with rotor and adapter for 50 mL centrifuge tube (non-skirted) •
- Microcentrifuge for 1.5 mL tube •
- Micropipettes (0.5-10 µL, 2-20 µL, 10-100 µL, 100-1000 µL)
- Freezer/Refrigerator (-20°C/4°C)
- Vortex Mixer
- Water Bath
- Additional Accessories d)
 - Tube rack/stand
 - Tweezers



Beveled pipette tip

² QIAGEN (Cat. No. 69514) ³ QIAGEN (Cat. No. 204741/204743)

⁴ National Scientific Supply Company, Inc. (Cat. No. BUN020GL-MRS) or VWR (Cat.No. 10126-388) is recommended or equivalent



Warnings and Precautions

- All samples should be treated as potentially infectious.
- Wear appropriate personal protective equipment, including (but not limited to) protective disposable gloves, laboratory coats and eye protection when handling samples and kit reagents. Wash hands thoroughly after handling specimens and reagents.
- Clean and decontaminate work area and instruments, including pipettes, with commercially available decontamination products.
- A designated working area should be dedicated for processing specimens and to add extracted specimens to PCR Reaction Mix.
- Use sterile pipette tips with filters.
- Handle Chip with care, avoid contact with the microreactor.
- Each Chip is used to process one test. Do not reuse processed Chips.
- Do not use kit, Chips or reagents after the expiration dates shown on the respective labels.
- Follow laboratory safety rules and procedures as defined by approved biohazard safety guidelines or regulations.
- Discard waste according to your local safety regulations.
- Material Safety Data Sheets (MSDS) are available upon request.

Quality Control

Under Veredus' quality assurance program, the performance of our VereFoodborne[™] Detection Kit is monitored routinely to ensure consistent product quality. Sampling is done on each lot and tests carried out via amplification of the respective control nucleic acid fragment.



Workflow





Specimen Collection, Handling and Storage

Specimens should be collected, handled and stored following the user institution's standard procedures. Inadequate or inappropriate specimen collection, storage and transport are likely to yield false negative results. Training in specimen collection is highly recommended because of the importance of specimen quality.

Store extracted nucleic acids at 2°C to 8°C and use within 24 hours. For longer storage, store extracted nucleic acids at -25°C to -15°C. Repeated thawing and freezing may affect the quality of the nucleic acid. They should be tested before use.

Sample Preparation

Extracted DNA is the starting material for the VereFoodborne[™] Detection Kit. To obtain maximum performance, it is very important to establish the extraction process. Some naturally-occurring substances, such as heme, melanin, and polysaccharides could be PCR inhibitors and could interfere with the assay performance. Please refer to the respective manufacturer's handbook for detailed extraction procedure.

For sample elution, EDTA-containing buffer (TE, with 0.1-1.0 mM EDTA) is the standard in most of the extraction kit. However, EDTA in the sample may inhibit the PCR process if it is used at higher (i.e. 10 mM) concentration. Please check the extraction kit components for EDTA concentration.

The following nucleic acid extraction kit is recommended:

• DNeasy *mericon* Food Kit⁵

⁵ QIAGEN (Cat. No. 69514)



Protocol

Important notes before starting

- Thaw all frozen components thoroughly at room temperature before use.
- After thawing, briefly mix and centrifuge the components.
- Vortex briefly or pipette up and down 5-6 times when mixing reagents with enzymes. Avoid making bubbles.
- Precipitation or crystallization may occur in Wash Buffer Concentrate if the storage temperature is low. Should this occur, please refer to important note before Chip washing step (page 23).
- User intervention is required:
 - After the PCR protocol is completed, user to load hybridization mix into the Chip and return Chip to the TCM for hybridization step.
 - After hybridization protocol is completed, user to wash the Chip and place the Chip into the OR for detection.
- Use current VereFoodborne[™] Chip version as follows:
 - FBN-30
- For software, use current version or higher as follows:
 - VerePLEX[™] Biosystem Version 5.3.X⁶
 - BioApplication
 VereFoodborne_3.0.7
- Screenshots are for illustration purposes only, and individual installations may vary.

⁶ Current VerePLEX[™] Biosystem version. Subject to minor changes.



1. Dilute sL PCR Control

i) Dilute the sL PCR Control (5x10⁴ copies/µL) to 5x10² copies/µL in a sterile PCR tube according to the following pipetting scheme:

Components	Volume (µL)
sL PCR Control	2
PCR Grade Water	198
Total	200

- ii) Mix thoroughly by pipetting or brief vortexing.
- iii) Briefly centrifuge to pull contents down to bottom of tube.

2. Prepare PCR Reaction Mix

i) Depending on the number of samples, prepare the required volume of PCR Reaction Mix in a sterile PCR tube according to the following pipetting scheme:

Number of Reactions	1 rea	1 reaction		
PCR Reaction Mix Component	Tube A (µL)	Tube B (µL)		
QuantiTect Multiplex PCR NoR	6.25	6.25		
	A	1	-	
	В	-	1	
Diluted sL PCR Control (5x10 ² o	copies/µL)	1	1	
Extracted Sample	2.5	2.5		
PCR Grade Water	1.75	1.75		
Total		12.5	12.5	

- ii) Mix thoroughly by pipetting or brief vortexing.
- iii) Briefly centrifuge to pull contents down to bottom of tube.



- 3. Load PCR Reaction Mix into Chip
- i) Insert a VereFoodborne[™] Chip into the Chip Holder and ensure a secure hold.
- NOTE: To ensure a secure and firm positioning of the Chip, the holder has a pin fastener, press to release holder (*Figure 1*) when inserting and removing the Chip.





ii) Draw **11.5 µL** of the PCR Reaction Mix **Tube A** with a pipette.

NOTE: Use a 20 μ L pipette and recommended pipette tips⁷ for Chip loading.

- iii) Hold the pipette in a vertical position, in such a way that the tip is perpendicular to the surface.
- iv) Fit the tip into one of the inlet holes (see Figure 2a).
- v) Applying slight pressure onto the tip, press the plunger smoothly to the first stop position (see *Figure 2b*), allowing the mix to flow into the PCR chamber.



Do not press the plunger <u>beyond the first stop</u> as this will introduce air into the chamber and mix will flow into the microarray chamber. Keep the plunger at the first stop until you remove the tip from the inlet (this procedure avoids spilling and the injection of air inside the chambers).

⁷ National Scientific Supply Company, Inc. (Cat. No. BUN020GL-MRS) or VWR (Cat.No. 10126-388) is recommended or equivalent



vi) Using a new pipette tip, repeat steps (ii) to (v) for another chamber. Load the Reaction Mix **Tube B** into the other inlet.



Figure 2: (a) Tip placement into an inlet hole during sample loading; (b) Pipette sketch indicating the different stop positions

Press and release the pipette plunger <u>slowly</u> at all the times. Never allow the push button to snap back. Check for foreign particles in the tip. Hold the pipette in an upright position while aspirating liquid.

4. Seal Chip for PCR

i) The IN and PCR sealing clamps are shown in *Figure 3*. The IN clamp (labeled "**2** IN") is dedicated to seal the inlet holes, and the PCR clamp (labeled "**1PCR**") to seal the outlet holes in the microarray chamber.



Figure 3: IN and PCR clamps

- ii) The clamp undersides are different, owing to their specific sealing function:
 - The "2 IN" clamp (*Figure 4*a) has an elastomer with a rectangular protrusion that seals the inlet holes, and one alignment pin that fits the corresponding hole on the Chip (*Figure 5*);



• The "**1PCR**" clamp (*Figure 4*b) has an elastomer with a rectangular protrusion that seals the outlet holes, and two alignment pins that fit the corresponding holes on the Chip (*Figure 5*).



Figure 4: Bottom view of (a) "2 IN" and (b) "1PCR" clamps



Figure 5: Alignment holes for "2 IN" and "1PCR" (or Hybridization) clamps

iii) Attach the "**1PCR**" first by pressing the lateral flyers and placing the pins into the alignment holes (*Figure 4* and *Figure 5*). After reaching the final position (the solid part of the clamp touches the edge of the Chip) release the flyers and press the upper part of the clamp until a 'click' sound is heard (*Figure 6*).



Figure 6: "1PCR" clamp attached to Chip



- iv) Repeat step (iii) with "2 IN" clamp.
- v) After the Chip is sealed (*Figure 7*), remove Chip from the Chip Holder.



Figure 7: Sealed Chip ready for PCR

5. Run Chip (PCR)

- i) Switch on the TCS.
- ii) Switch on the computer and launch "E@syControl" software by clicking on the icon

son the computer desktop.

The program will start searching for connected TCSs and the green TCS icon will be displayed in the "*E@syControl*" window when TCS is connected. TCM will display "*READY TO USE*" message on the LCD screen.

iii) Click "*Login*" on the toolbar. The "*Login*" window will be displayed. Log in with the correct username and password.

E@syCon	rol	Help			
Clogin	Configure BioApplicatio	on 🖉 Platform 🛃 B	ixit		
				en e	@syControl
TC and	S 236				۲
тсм	Barcode/Chip Id	BioApplication	State	Actip	User
🥟 ТСМ 1	22123E0H06012	Default_BioApplica v	0%	togin St	ſ
🥟 ТСМ 2		Select a program v	0%	s E@syControl	ſ
🥟 ТСМ 3		Select a program v	0%	St Insert usemame and password	ſ
🥟 ТСМ 4		Select a program v	0%	St Username: user	ſ
🥟 ТСМ 5		Select a program 👻	0%	St Password:	La C
				Login Cancel	



iv) Select "*BioApplication*" from the toolbar and check the "**VereFoodborne**" (version 3.0.7 or higher) BioApplication.

figuratio	on To onfigure	ols Platform Help	Platform					
Bi	oApplicat	tion Management			X		E@syCo	nt
Ba	Default	BioApplication Name	Temperature Program	Created by			User	
		Default_BioApplication VereFever 1.0.1	Prot_default_1 FVR_UDG_RT_20	sysadmin sysadmin	_	ip insertion	easypowe	r i
		VereFlu_3.8.1 VereFoodborne_3.0.6	Fu_3.8.1	sysadmin		and one set of calibration data	easypowe	r
j		VereMERS_1.1.1	MERS_MT	sysadmin		and one set of calibration data	easypowe	r u
		VereThreat_2.4.1	Threat_2.4.1	sysadmin		and one set of calibration data	easypowe	r i
		vere vet_1.3.2	vere vet_1.3.2	sysadmin		and one set of calibration data	easypowe	r l
					Ok			
L					Cancel			

- v) Open the lid of the TCM (if not already open).
- vi) Place Chip into TCM.
- NOTE: Ensure the alignment pins on the TCM are inserted into the corresponding alignment holes on the Chip (*Figure 8*).



Figure 8: Chip inserted into the TCM



vii) In the "*E*@syControl" window, select the appropriate TCS and place the cursor in the "*Barcode/Chip Id*" field of the respective TCM.

👒 E@syCont	rol	_	-			- 0 -	x		
File Con	figuration Tools Platform	Help							
🔆 😡 Logout 🙊 Configure 🗖 BioApplication 🛷 Platform 📓 Exit									
					en 💭 👘	@syCon	trol		
E IC	S 236						۲		
тсм	Barcode/Chip Id	BioApplication	State	Action	Status	User			
🛷 ТСМ 1	22123E0H06012	Default_BioApplica -	0%	Start	Not Ready - Wait for correct chip insertion	easypower	a		
🛷 ТСМ 2		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ		
🧳 ТСМ З		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	<u> </u>		
🛷 ТСМ 4		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ		
🧳 ТСМ 5		Select a program 👻	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	<u>_</u>		

- viii) Scan the 2-D barcode on the Chip using the barcode scanner.
- NOTE: Wait for the program to register the Chip calibration data before scanning the next Chip. The Chip ID should become "green" color.

👒 E@syCont	trol			-		- 0	x	
File Con	figuration Tools Platform	Help						
😔 Logout 🙊 Configure 🗖 BioApplication 🎻 Platform 👼 Exit								
🚎 E@syCa								
E TC	۲CS 236							
TCM	Barcode/Chip Id	BioApplication	State	Action	Status	User		
🧳 ТСМ 1	22123E0H06012	Default_BioApplicar 💌	0%	Start	Not Ready - Wait for correct chip insertion	easypower	<u> </u>	
🔗 ТСМ 2	2210110R20038	Default_BioApplicar •	0%	Start	Not Ready - Wait for correct chip insertion	easypower		
🧳 ТСМ З		Select a program 💌	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ	
🛷 ТСМ 4		Select a program 👻	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ	
🛷 тсм 5		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ	

- ix) Repeat step (v) to (viii) to register all of the Chips to be run.
- x) For random access, select the relevant BioApplication from the dropdown list under *"BioApplication"* for each TCM.

File Cont	figuration Tools Platf	form Help						
O Logout	Configure 🔝 Option	ns 🛛 🥅 Temperature Program	BioAppl	ication 🖄 G	raphs 👻 🖉 Platform 🛛 🌉 Exit			
				- Hereit		æ .	@sy Co	ntro
E ICS	5 078							(
TCS	6 098							(
тсм	Barcode/Chip Id	BioApplication	State	Action	Status		User	
🖉 ТСМ 1	22128H0A13130	Default_BioApplica 👻	0%	Start	Not Ready - Wait for correct chip insertion		admin	4
🖉 ТСМ 2	22127GDM16006	Default_BioApplica 👻	0%	Start	Not Ready - Wait for correct chip insertion		admin	1
🖉 ТСМ З		Select a program Default_BioApplication	0%	Start	Not Ready - Wait for program and one set of calibration data		admin	1
🖉 ТСМ 4		TestOPQ1 TestOPQ2	0%	Start	Not Ready - Wait for program and one set of calibration data		admin	6
A TCM 5		VereFlu_3.3	0%	Start	Not Ready - Wait for program and one set of calibration data		admin	3



xi) After selection, right-click the selected BioApplication. An option menu will be displayed. Select "Send" to load the selected program to the corresponding TCM.

🛪 E@syCont	trol		_			- 0 -	x
File Con	figuration Tools Platform	Help					
🥥 Logout	🛛 👰 Configure 🛛 🗖 BioApplicati	ion 🛷 Platform	🧕 Exit				
					<u>بین</u> E	@syCor	itrol
E IC	S 236						۲
тсм	Barcode/Chip Id	BioApplication	State Action	n	Status	User	
🥔 ТСМ 1	22123E0H06012	Default_BioApplica	• 0% Start		Not Ready - Wait for correct chip insertion	easypower	a
🧳 ТСМ 2	2210110R20038	Default_BioApplica	Start		Not Ready - Wait for correct chip insertion	easypower	ſ
🧳 ТСМ 3		Select a program	Save temp profile		Not Ready - Wait for program and one set of calibration data	easypower	ſ
🛷 ТСМ 4		Select a program	Send		Not Ready - Wait for program and one set of calibration data	easypower	ſ
🥔 ТСМ 5		Select a program	Send group View		Not Ready - Wait for program and one set of calibration data	easypower	ſ
			Security Lock Lock all Advanced View log Reset				

- xii) Close the lid of TCM. The TCM will validate the Chip against its calibration data and the TCM will display "CHIP VALIDATION" message on LCD screen.
- xiii) Once the Chip is validated, the TCM will display "CHIP INSIDE PRESS START" on the LCD screen or "Ready" in the "Status" field in the "E@syControl" window.
- xiv) Press "*Play*" button (▶) on the TCM front panel or click "*Start*" in the "*E*@syControl" window to begin thermal program.

👒 E@syCont	rol	_		-			x		
File Con	figuration Tools Platform	Help							
😔 Lagout 🙊 Configure 🗖 BioApplication 🛷 Platform 🧕 Exit									
🚎 E@syCo									
E IC	TCS 236								
TCM	Barcode/Chip Id	BioApplication	State	Action	Status	User			
🛷 ТСМ 1	22123E0H06012	Default_BioApplica -	0%	Start	Not Ready - Wait for correct chip insertion	easypower	ſ		
🔗 ТСМ 2	2210110R20038	Default_BioApplica -	0%	Start	Ready	easypower	ſ		
🛷 ТСМ З		Select a program 👻	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ		
🥟 ТСМ 4		Select a program 💌	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ		
🛷 ТСМ 5		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ		

xv) Once the PCR protocol is completed, "WAITING FLUIDIC OPERATION" will be displayed on the LCD screen of the TCM or "Fluidic Operation – Wait for user to open the lid" will be displayed under the "Status" field in the "E@syControl" window.

S E@syControl File Configuration Tools Platform Help										
🥥 Logout	© Logout ☆ Configure ☐ BioApplication									
E TC	S 236					e ,	٢			
TCM	Barcode/Chip Id	BioApplication	State	Action	Status	User				
🛷 ТСМ 1	22123E0H06012	Default_BioApplicar 💌	0%	Start	Not Ready - Wait for correct chip insertion	easypower	ſ			
🔗 ТСМ 2	2210110R20038	Default_BioApplica 👻	87%	Stop	Fluidic Operation - Wait for user to open the lid	easypower	_			
🛷 тсм з		Select a program 👻	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	<u>_</u>			
🛷 ТСМ 4		Select a program 💌	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	_			
🛷 тсм 5		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	£			



6. Prepare Hybridization Mix

i) Prepare Microarray Hybridization Mix by transferring **870 μL** of VCP Hyb Buffer to 1 tube of **30 μL** Hyb Probe Concentrate:

Number of Reactions	25 reactions
Hyb Mix Components	Volume (µL)
VCP Hyb Buffer	870
Hyb Probe Concentrate	30
Total	900

ii) Leave Microarray Hybridization Mix at room temperature to equilibrate for at least 20 minutes.



Microarray Hybridization Mix must be equilibrated at room temperature for 20 mins before use. Mix well before use.

iii) Mix thoroughly by vortexing the tube briefly (~10 seconds) or inverting it 4-6 times before spinning down.

7. Load Hybridization Mix into Chip

- i) Remove Chip from the TCM when prompted.
- ii) Insert Chip onto the Chip Holder and ensure a secure fit.
- iii) Remove "2 IN" and "1PCR" clamps and discard them. DO NOT reuse the clamps.
- iv) Draw **14.5 µL** of the Microarray Hybridization Mix with a pipette.

NOTE: Use a 20 µL pipette and recommended pipette tips⁸ for Chip loading.

- v) Hold the pipette in a vertical position, in such a way that the tip is perpendicular to the surface.
- vi) Fit the tip into one of the inlet holes (see *Figure 9a*).

⁸ National Scientific Supply Company, Inc. (Cat. No. BUN020GL-MRS) or VWR (Cat.No. 10126-388) is recommended or equivalent



Applying slight pressure onto the tip, press the plunger smoothly to the first stop vii) position (see Figure 9b), allowing the mix to flow into the PCR chamber. The PCR mix inside the PCR chamber will be displaced by the Microarray Hybridization Mix and will be observed to fill up the microarray chamber (Figure 10).



Do not press the plunger beyond the first stop as this will introduce air into the chamber and mix will flow into the microarray chamber. Keep the plunger at the first stop until you remove the tip from the inlet (this procedure avoids spilling and the injection of air inside the chambers).

viii) Using a new pipette tip, repeat steps (vi) and (vii) for another chamber. Load the mixture into the other inlet.

NOTE: Use a new tip for every loading to prevent carryover of the PCR product.



Figure 9: (a) Tip placement into an inlet hole during sample loading; (b) Pipette sketch indicating the different stop positions

(b)



Press and release the pipette plunger slowly at all the times. Never allow the push button to snap back. Check for foreign particles in the tip. Hold the pipette in an upright position while aspirating liquid.



ix) Tap the Chip gently at the side if the solution in the microarray chamber does not fully fill the microarray chamber.



Figure 10: *Microarray chamber filling* – (*a*): *filling of first inlet;* (*b*) *filling of second inlet* (*c*): *completely filled*

- 8. Seal Chip for Hybridization
- i) Prepare a new IN clamp (labeled "2 IN") and HYB clamp (not labeled) as shown in *Figure 11*.



Figure 11: IN and HYB clamps

ii) The HYB clamp has a flat PDMS surface, 100 µm deep elicited in the gasket, and a surrounding trench used to accommodate the air displaced by the solution (*Figure 12*).



Figure 12: Bottom view of HYB clamp



iii) Attach the "2 IN" clamp first (*Figure 13*).



Figure 13: "2 IN" clamp attached to Chip

iv) Seal the microarray chamber carefully using the HYB clamp, making sure that no bubbles are introduced into the chamber (*Figure 14*).



Figure 14: Sealed Chip ready for the hybridization

NOTE: Should any bubbles form during the sealing of the microarray chamber, tap the Chip gently on the workbench, microarray closest to the bench surface. This will force the bubbles to migrate to the outlet edge of the PCR chamber where there are no probes (*Figure 15*).



Figure 15: Outlet edge of the PCR chamber on the Chip



v) Remove Chip from the Chip Holder, making sure the clamps are tightly held in place.

9. Run Chip (Hybridization)

i) Load the sealed Chip into the respective TCM. Press "*Play*" button (▶) on the TCM front panel or press '*Start*' when prompted.

🐞 E@syCont	trol		-	-			x	
File Configuration Tools Platform Help								
😔 Logout 🙊 Configure 🗖 BioApplication 🔗 Platform 🛃 Exit								
🚎 E@syCo								
7CS 236								
тсм	Barcode/Chip Id	BioApplication	State	Action	Status	User		
🥟 ТСМ 1	22123E0H06012	Default_BioApplica 👻	0%	Start	Not Ready - Wait for correct chip insertion	easypower	a	
🛷 ТСМ 2	2210110R20038	Default_BioApplica 👻	87%	Start	Fluidic Operation - Sensor valid - waiting for a start	easypower	_	
🧳 ТСМ 3		Select a program 👻	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ	
🔗 ТСМ 4		Select a program 💌	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ	
🛷 ТСМ 5		Select a program 👻	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ	

ii) Once the hybridization protocol is completed, "COMPLETED" message will appear on the LCD screen of the TCM or "Completed – Wait for user to open the lid" will be displayed in the "E@syControl" window. Remove the Chip from TCM and proceed to washing step (Section 10) immediately.

r& E@syControl File Configuration Tools Platform Help ♀ Logout I <a>Configure ■ BioApplication ♪ Platform ■ Loit										
	E@syControl									
TC:	5 236						۲			
ТСМ	Barcode/Chip Id	BioApplication	State	Action	Status	User				
🥟 ТСМ 1	22123E0H06012	Default_BioApplicar •	0%	Start	Not Ready - Wait for correct chip insertion	easypower	ſ			
🧳 ТСМ 2	2210110R20038	Default_BioApplica -	100%	Start	Completed - Wait for user to open the lid	easypower	ſ			
🛷 тсм з		Select a program 👻	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ			
🛷 ТСМ 4		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	ſ			
🛷 тсм 5		Select a program 🔹	0%	Start	Not Ready - Wait for program and one set of calibration data	easypower	a			



Temperature is one of the biggest factors to control hybridization process. TCM will not maintain the temperature after hybridization process is completed. To avoid exposure to the lower temperature, it is highly recommended to start washing step immediately after the hybridization protocol is completed.

10. Wash Chip

In case of precipitation or crystallization in the Wash Buffer Concentrate, warm up the entire bottle of buffer in a water bath set to 42°C for at least 1 hour with occasional shaking until all the precipitate is dissolved. Mix thoroughly before use.



i) Measure 50 mL of the supplied Wash Buffer Concentrate into a 1 L empty bottle. Top up the 1 L bottle with Distilled/Reverse Osmosis (RO)/Ultrapure Water to 1 L. Mix well.

Components			Volume (mL)
Wash Buffer Conce	50		
Distilled/Reverse	Osmosis	(RO)/Ultrapure	
Water			950
Total			1000

ii) Prepare and fill non-skirted 50 mL centrifuge tubes with 50 mL of the prepared Wash Buffer from step (i) (Figure 16).

In case of precipitation or crystallization in the diluted Wash Buffer, warm up the entire bottle of buffer in a water bath set to 42°C with occasional shaking until all the precipitate is dissolved. Equilibrate at room temperature for at least an hour and mix thoroughly before use.

- NOTE: Fill Wash Buffer to 50 mL mark on centrifuge tube to completely submerge the Chip. Place only ONE Chip per tube.
 - Remove the "2 IN" and HYB clamps, paying attention that no liquid spills out. Discard iii) clamps. **DO NOT** reuse the clamps.
 - iv) Insert Chip with the microarray end at the top into the centrifuge tube (Figure 17). Screw the tube cap on.







Figure 16

Figure 18

V) Place the tube with the microarray side facing towards the rotor axis (Figure 18).



vi) Centrifuge the tube at 3000 rpm for 2 minutes.



The centrifuge spins at high speeds. Ensure that the lid is closed properly and that all the buckets are correctly balanced.

- vii) After centrifugation, empty the tube of Wash Buffer and place tube back into the centrifuge with the microarray in the same orientation as step (v).
- viii) Centrifuge the tube at 3000 rpm for 2 minutes to spin-dry the microarray.
- ix) After centrifugation, remove the Chip using a pair of tweezers.
- x) Proceed to detection step immediately.



Fluorescent dye is used and is prone to degradation upon ozone exposure. It is highly recommended to proceed to the detection step immediately after washing. Minimize exposure of hybridized arrays to light, high temperatures and high ozone levels after washing. Place the Chip with the microarray face down on a clean paper towel or in a container.

11. Detection

- i) Switch on the Optical Reader (OR).
- ii) Launch "*E@syCheck*" software by clicking on the icon whe computer desktop.
- iii) Click "*Login*" on the toolbar. The "*Login*" window will be displayed. Log in with the correct username and password.

ports						
BarCode	Report Name	BioApplication	Status	User	Note	State
		(Lucia	×			
		Login				
			🔮 E@syCheck			
		Insert username and pass	E@syCheck			



iv) Click "BioApplications" on the toolbar.

& E@syCheck - [New Report]	E@syCheck - [New Report] By File View Tools Platform Window Help										
🙆 Logout 🔄 New 💕 Open	🕒 Logout 🐘 Hevr 🔊 Open 🖫 Save 🐇 Print 🛝 Acquisition 😿 BoAgokaston 🤣 Platform 🛃 Ext										
Reports											
BarCode	Report Name	BioApplication	Status	User	Note	State					

v) The "*BioApplication Management*" window will be displayed. Check the "**VereFoodborne**" (version 3.0.7 or higher) BioApplication. Click "*Ok*" to proceed.

Name	User	Default	
Default_BioApplication	sysadmin	 Image: A start of the start of	
VereFever_1.0.1	sysadmin		
VereFlu_3.8.1	sysadmin		
VereFoodborne_3.0.6	sysadmin		
VereMERS_1.1.1	sysadmin		
VereMTB_4.3.1	sysadmin	1	
VereThreat_2.4.1	sysadmin		
VereVet_1.3.2	sysadmin		

vi) Click on the "Barcode" field to bring the cursor to this location.

& E@syCheck - [New Report]										
🚯 File View Tools Platform Window Help										
🕒 Logout 🖹 New 🗊 Open 🔄 Save 🍊 Frint 🔱 Acquisition 💥 BioApplication 🔗 Platform 🛃 Exit										
Reports										
BarCode	Report Name	BioApplication	Status	User	Note	State				

vii) Scan the 2-D barcode on the respective Chip.



viii) Open the OR lid (if not already open).



ix) Insert Chip into the OR with the microarray facing up (*Figure 19*).



Figure 19: Chip inserted into OR

x) Click on "*Acquisition*" button on the toolbar. The OR will begin image acquisition and image analysis immediately.

& E@syCheck - [New Report]					-					
🕼 File View Tools Platform Window Help										
🕝 Logout 🐘 New 💰 Open 🔄 Save 🍊 Frint 🔖 Acquisition 🛠 BoAppleation 🛷 Platform 🛃 Ext										
Reports										
BarCode	Report Name	BioApplication	Status	User	Note	State				
2210110R20038		Default_BioApplication		easypower						

xi) After "*Acquisition*" operation, "*Completed*" will be displayed in the "*Status*" field and the results will be displayed.





Clicking the "Final Check" tab shows target/control detection summary.

t Result Final Check	
Target	Result
T_100	Not Detected
T_200	Detected
T_300	Not Detected
T_400	Not Detected
T_500	Not Detected
T_600	Not Detected
T_700	Not Detected
T_300 / T_400	Not Detected
T_800 / T_900	Not Detected

xii) Comments on this particular Chip run can be recorded under the "*Note*" field and this will be printed on the final report.

\$	E@syCheck - [New Report]										
	File View Tools Platform	Window Help						_ 8 ×			
0	🕒 Logout 🐘 New 🗊 Open 📓 Save 👼 Print 🚯 Acquisiton 💥 BoApplication 🛷 Platform										
	Reports										
	BarCode	Report Name	BioApplication		Status	User	Note	State			
	2210110R20038		Default_BioApplication	Completed	66	sypower					

- NOTE: "Note" field can be modified only BEFORE saving the analysis.
 - xiii) To save the analysis, press "*Save*" button on the toolbar. All the information associated to the analysis will be stored in a local database.

🗞 E	🗞 E@syCheck - [New Report]											
В.	🚯 File View Tools Platform Window Help											
۵	🕒 Logout 🖹 New 🗊 Open 📆 Save 👼 Print 🗞 Acquisiton 💥 BioApplication 🧳 Platform 📲 Ext											
	Reports											
	BarCode	Report Name	BioApplication		Status	User	Note	State				
	2210110R20038		Default_BioApplication	Completed		easypower						

xiv) Click "File" on the top toolbar and select "Print Report" from the context menu.

\$	E@sy	Check - [New Report]				Contraction from the		
	File	View Tools Platform	Window Help					_ 8 ×
6		New Report	Save 进 Print 🔹 Acquisition 💥 BioAp	oplication 🔗 Platform 🚮 Exit				
		Open Report						
		Save Report	Report Name	BioApplication	Status	User	Note	State
		Print Report						
	\checkmark	BioApplication	0110R-20-38	Default_BioApplication	Completed	easypower		
		Exit						



xv) The "*Print Report*" window will be displayed. Click "*Print*" to proceed.

NOTE: By default, the microarray equivalent image is not included in the PDF report, check "*Print Microarray*" to print report with the microarray equivalent image.

Options Print Microarray Customer Logo	Page Settings Portrait
File	Browse

xvi) A "Save PDF file as" window will be displayed. Save file in the desired destination folder. Click "Save" to proceed.



xvii) Alternatively, press "*Print*" button on the toolbar to display and print the PDF report.

🗞 E(DsyCheck - [New Report]						
ا 🛃	🔉 File View Tools Platform Window Help						
O L	ogout 🛛 🗟 New 💕 Open	📲 Save <i>曇</i> Print 🛛 🗞 Acquisition 💢 BioApp	olication 🔗 Platform 🛃 Exit				
F	Reports						
	BarCode	Report Name	BioApplication	Status	s User	Note	State
	2210110R20038		Default_BioApplication	Completed	easypower		



Assay Control

The following controls are included in each test:

- 1. Internal PCR controls to check for succesful nucleic acid amplification reaction
- 2. Positive and negative hybridization controls to check for hybridization-related issues
- 3. Orientation controls for microarray grid alignment

It is recommended to include a negative control sample for each test run to check for possible contamination.

Interpretation of Results

The VerePLEXTM Biosystem software provides a qualitative result for the presence (Detected) or absence (Not Detected) of the target gene/organism under the "*Test Result*" and "*Final Check*" tab on the software interface and the final printed report. "Inconclusive" will be displayed for presence of bad spot signal(s) and the number of bad spot fails to meet the set criteria.

The control fields will display "Valid", "Not Valid" and "Inconclusive". If the hybridization control and/or and negative control is "Not Valid" or "Inconclusive", the software will not proceed with any further data analysis and the result of the control will only be displayed in the "*Test Result*" section. No result will be displayed in the "*Final Check*" section. The validity of PCR control has no influence of the outcome of the result if the target is "Detected".

For further details regarding the interpretation of the results and recommended actions, please refer to Troubleshooting Guide or contact our Technical Support.



Troubleshooting Guide

The troubleshooting guide may be helpful in solving problems that may arise.

Comments and Recommended Actions		
Instrument and Software Issues		
Make sure that the instruments (TCS and OR) are properly maintained. Refer to the Troubleshooting section of the VerePLEX [™] Biosystem SYSTEM & SOFTWARE IFU.		
Refer to the Troubleshooting section of the VerePLEX™		
Biosystem SYSTEM & SOFTWARE IFU.		
 Make sure the correct BioApplication is used. One or more negative control probes have a fluorescent signal due to a high fluorescent background or fluorescent artifacts in the microarray surface. This problem may occur in the following situations: Dust or fiber-like material is found on the microarray area. Inspect the area and repeat detection step. If the high background persists repeat the test. The Chip is not washed properly. Repeat the washing step following carefully the Instructions for Use (IFU). If the high background persists repeat the test. The Wash Buffer is not diluted according to the instructions. Prepare a new bottle of Wash Buffer following carefully the IFU and repeat the washing step using the newly prepared Wash Buffer. If the high background persists repeat the test. The HYB clamps are not properly inserted causing partial or total evaporation of the Microarray Hybridization Mix with a consequent drying of the fluorescent mix on the microarray surface. Repeat the test. Contamination. Repeat the test. Ensure that the workspace and instruments are decontaminated at regular intervals. Refer to the Cleaning section of the respective IFU. 		
1. Make sure the correct BioApplication is used.		



	Comments and Recommended Actions
valid" OR "Inconclusive"	 Some of the hybridization control probes are not lit up. This problem may occur when the labeled hybridization probes are partially degraded due to wrong storage conditions. Repeat the test using new reagents and follow carefully the IFU on Storage Condition. High fluorescent background or fluorescent artifacts in the microarray surface (see negative controls not valid or inconclusive comments and suggestions).
3. No signal	 Make sure the correct BioApplication is used. Photo-bleaching of the dye signal in the microarray detection area due to high level of ozone in the lab. Repeat the test and minimize the microarray surface exposure to the light as much as possible. The labeled primers and hybridization probes are degraded due to wrong storage conditions. Repeat the test using new reagents and follow carefully the IFU on Storage Condition. The Wash Buffer is not diluted according to the instructions. Prepare a new bottle of Wash Buffer following carefully the IFU and repeat the washing step using the newly prepared Wash Buffer. Defective Chip. Repeat test with new Chip. The kit has expired. Check the expiry date of the kit and use a new kit, if necessary.
PCR control "Not Valid"	
 No fluorescent signals for PCR control probes but the specific target probes give signals 	 Make sure the correct BioApplication is used. This may occur when the target nucleic acid is much more concentrated than the PCR control. This has no influence on the outcome of the test. The PCR control is degraded due to wrong storage conditions. This has no influence on the outcome of the test but use new reagents in the next run.
 No fluorescent signals for both the PCR control probes and the specific target probes 	 Make sure the correct BioApplication is used. The PCR control and the target nucleic acid are degraded due to wrong storage conditions. Repeat the test using new reagents and sample and follow carefully the IFU on Storage Condition. The PCR control and the target nucleic acid are degraded due the presence of DNase. Check quality of nucleic acid sample or use fresh nucleic acid sample.



	Comments and Recommended Actions			
	 Repeat the test and follow carefully the IFU on Warnings and Precautions. 4. The labeled primers are degraded due to wrong storage conditions. Repeat the test using new reagents and follow carefully the IFU on Storage Condition. 5. PCR was inhibited. Use recommended extraction kit. Refer to the manufacturer's handbook for detailed extraction procedure. Repeat the test. 6. Defective Chip. Repeat test with new Chip. 7. The kit has expired. Check the expiry date of the kit and use a new kit, if necessary. 			
Inconclusive results for specific probes				
 One or more specific probes have bad spots 	Some of the spots of one or more specific probes are not recognized as a spot by the E@syCheck software. This problem may occur when the spot morphology is not good mainly due to a high fluorescent background or fluorescent artifacts on the microarray surface (see negative controls not valid or inconclusive comments and suggestions).			
2. One or more probes have one replica of a specific spots pair not lit up	 This problem may occur when there is not enough amplified dye-labeled target to hybridized the microarray due to the following reasons: 1. The labeled primers are partially degraded due to wrong storage conditions. Repeat the test using new reagents and follow carefully the IFU on Storage Condition. 2. The target nucleic acid is partially degraded due the poor sample preparation and storage. Check quality of nucleic acid sample or use fresh nucleic acid sample. Repeat the test and follow carefully the IFU on Warnings and Precautions. 3. Insufficient starting material. Repeat test with increased amount of nucleic acid sample. 			



Limitations of the Test

- Use of this kit should be limited only to trained personnel.
- This test is a qualitative test and does not provide a quantitative value for the detected pathogen in the sample.
- Strict compliance with the IFU is required for optimal results. Modifications to these procedures may alter performance of the test.
- Appropriate specimen collection, handling, storage and processing procedures are required for the optimal performance of this test.
- This test is not to be used on specimen directly. Specimen needs to be processed using appropriate nucleic acid extraction methods prior to using this test.
- The dye used is susceptible to degradation upon exposure to ozone. Strict compliance with the processing procedures is required for optimal performance of the test. If possible, procedures should be done in a reduced ozone environment to eliminate degradation of the dye molecule.
- It is advised to scan each microarray only once, subsequent scans may not yield similar results as fluorescence intensity may decrease due to decay of the fluorophore.
- Results from the test should be interpreted with other laboratory data available.
- Although the kit is highly specific and sensitive, a low incidence of false results can occur. A negative result does not preclude the possibility of existence of the target organisms in the sample. Other available tests are required if questionable results are obtained.
- A specimen yielding a negative result may contain other organisms other than the target organisms.
- False negative results may occur due to presence of sequence variants in the gene targets of the assay, procedural errors, amplification inhibitors in specimens, or inadequate nucleic acids for amplification.
- False positive results may occur due to cross-contamination by target organisms, their nucleic acids, amplicons, or from non-specific signals in the test.
- Cross-reactivity with organisms not tested can lead to erroneous results.
- Bacteria nucleic acids may persist *in vivo* independent of organism viability. Detection of analyte target(s) do not imply that the corresponding organisms are infectious.



Disposal

Dispose of hazardous or biologically contaminated materials according to local safety regulations.

Technical Assistance

If you have any questions or technical issues regarding the use of the kit, or any other Veredus products, please contact our technical support department.

Contact

Your opinions, comments, questions or feedback are important to us and all Veredus' customers. Please contact us if you have any suggestions about product performance or new applications and techniques.

For information and technical assistance, please contact us via:

Telephone: +65 6496 8600

Fax: +65 6779 2680

Email: <u>info@vereduslabs.com</u>

Visit our website: www.vereduslabs.com



Understanding the Symbols

Symbol	Meaning
REF	Catalog number
LOT	Lot number
Σ	Contains sufficient for <n> tests</n>
2	Do not re-use
\sim	Date of manufacture
	Manufacturer
	Temperature limitation
\square	Use-by date (YYYY-MM-DD)
ī	Consult Instructions for Use
CONTROL +	Positive Control
Ţ	Fragile, handle with care
×	Keep away from sunlight
Ĵ	Keep dry
<u> </u>	This side up
Λ	Caution



Product Use Limitations, Warranty Disclaimer

Veredus Laboratories Private Limited (Veredus) manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products. Veredus products contain chemicals, which may be harmful if misused. Due care should be exercised with all Veredus products.

Veredus warrants that its products will meet the specifications stated on each product's specification sheet. If any component of the product does not conform to these specifications, Veredus will, at its sole discretion, as its sole and exclusive liability and as the users' sole and exclusive remedy, replace the product free of charge.

THIS WARRANTY LIMITS VEREDUS' LIABILITY TO THE REPLACEMENT OF THIS PRODUCT OR REFUND OF THE COST OF THE PRODUCT. NO OTHER WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT, ARE PROVIDED BY VEREDUS. VEREDUS SHALL HAVE NO LIABILITY FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL OR INCIDENTAL DAMAGES ARISING OUT OF THE USE, THE RESULTS OF USE OR THE INABILITY TO USE THIS PRODUCT AND ITS COMPONENTS.

In no event shall Veredus be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or use or the failure of Veredus products to perform in accordance with the stated specifications.

Some components of nucleic acid analysis, such as specific methods and compositions for manipulating or visualizing nucleic acids for analysis, may be covered by one or more patents owned by other parties. Similarly, nucleic acids containing specific nucleotides sequences may be patented. Making, using, offering for sale, or selling such components or nucleic acids may require one or more licenses. Nothing in this document should be construed as an authorization or implied license to make, use or sell any so covered component or nucleic acid under any such patents.



Notice to Purchaser

The procedure outlined in this protocol contains proprietary information. By purchasing this product, the user is granted a limited license by Veredus Laboratories Pte Ltd to use this information as described. The user of this product agrees not to use any of these proprietary methods in any other application and agrees not to communicate (either orally or in writing) these proprietary methods to any other person or institution.

VereFoodborne[™], VerePLEX[™] and VereChip[™] are trademarks of Veredus Laboratories Pte Ltd (Singapore, SG)

All other names of products and brands stated in this IFU are trademarks or registered trademarks of their respective owners.

© 2018 Veredus Laboratories Pte Ltd. All Rights Reserved.